

Sandata Mobile Connect

iOS and Android Global Release


Date: December 12, 2024

Faride Beaubien, Product Management



VC-1461: SMC ADA - Focus should not be able to move behind modals

- Both iOS and Android
- ISSUE: When the user executes any action that spawns a popover, drawer or modal, focus should not be able to reach items that are behind it.
 - Start unknown visit modal
 - Select a service popover
 - Select a location popover
 - o Visit Summary popover
 - Recipient Verification modal
- FIX: Focus will only select what is visually available to user






VC-1698: Front end pop to prevent users from logging in

frequently

ISSUE

- Both iOS and Android
- If SMC user attempts to login multiple times without logging out properly, the user will receive an unclear message
- Unclear message caused more failed attempts at logging in






VC-1698: Front end pop to prevent users from logging in

frequently

FIX

- Both iOS and Android
- If an SMC user attempts to login multiple times without logging out properly, the user will receive a clearer message
- Message: "You recently logged in on another device. Please try again after 2 mins."






VC-1710: SMC capturing GPS coordinates

ISSUE: Do not capture GPS coordinates if user changes device or app level permissions manually and logs into the app without quitting in the background

- Android and iOS
- Previous requirement to cache GPS coordinates in event there was a loss of signal

FIX:

- Location Services shall only be managed by the phone or SMC app level settings.
- Device settings will supersede SMC app settings. Location Services must be ON for any GPS coordinates to be captured AND SMC location permission must be granted.
- If user changes settings amidst a visit, do not capture GPS settings and ensure changes to settings are honored.









Scenario: Phone Location Services OFF; SMC Location Permission Changed ON -->OFF

- Turn on Location Services for device and allow location access to the SMC app
- Login and start any type of visit
- Go to device settings and change app location permissions to "Never" Do not quit app
- Login and finish the visit









- Both iOS and Android
- ISSUE: Existing labels for Chinese are mislabeled with the formal name instead of the written name
- FIX:
 - 。 Chinese (Mandarin) 普通話 -> Simplified Chinese (translation pending)
 - o Chinese (Cantonese) 粵語 -> Traditional Chinese (translation pending)







Sandata EVV

8.39 Release

Date: December 11, 2024

Sean Dunham, Product Management



HOM-7665: Add Payer, Program and Relationship to XREF Search Grid

- ISSUE: For Live-In Caregivers, there is a need to identify the payer, program and relationship when setting up XREF
 - Payer
 - Program
 - Live-In Caregiver Relationship
- FIX: The added columns will be displayed in the search results grid for the XREF module

